

ZNF433 Polyclonal Antibody



Catalog Number: 162379

Product name

ZNF433 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZNF433

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 564-673 of human ZNF433 (NP_001073880.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - ZNF433 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using ZNF433 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ZNF433 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using ZNF433 antibody at dilution of 1:100 (40x lens).





Immunohistochemistry - ZNF433 Polyclonal Antibody

ZNF433 Polyclonal Antibody

Catalog Number: 162379


Immunohistochemistry of paraffin-embedded human appendix using ZNF433 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - ZNF433 Polyclonal Antibody

Immunohistochemistry - ZNF433 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using ZNF433 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse lung using ZNF433 antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - ZNF433 Polyclonal Antibody

Immunohistochemistry - ZNF433 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using ZNF433 antibody at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded mouse brain using ZNF433 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZNF433 Polyclonal Antibody

ZNF433 Polyclonal Antibody



Catalog Number: 162379

Immunohistochemistry of paraffin-embedded

mouse kidney using ZNF433 antibody at dilution of 1:100 (40x lens).