

# MAGEA11 antibody



Catalog Number: 112379

## Product name

MAGEA11 antibody

## Specificity

Human; other species not tested.

## Antibody description

MAGEA11 Rabbit Polyclonal antibody. Positive IF detected in A375 cells. Positive WB detected in A375 cells. Observed molecular weight by Western-blot: 48 kDa, 44 kDa

## Preparation

This antibody was obtained by immunization of MAGEA11 recombinant protein (Accession Number: XM\_017029521). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:112379(MAGEA11 antibody) at dilution of



Immunofluorescent analysis of A375 cells using Catalog No:112379(MAGEA11 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)