

Anti-SPARCL1 / SPARC-like 1 antibody



Catalog Number: 104379

Product name

Anti-SPARCL1 / SPARC-like 1 antibody

Specificity

Human SPARCL1

Antibody description

Mouse monoclonal to SPARCL1 / SPARC-like 1

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SPARCL1 (rh SPARCL1; NP_004675.3; Met1-Phe664) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

FCM

Dilutions

FCM

Validations



SPARCL1 Antibody (FITC), Mouse MAb, Flow cytometric

Profile of anti-SPARCL1 reactivity on HL60 cells analyzed by flow cytometry. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat.No.554714), and stained with FITC conjugated Mouse anti-SPARCL1, The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.