

Anti-ZNF106 antibody



Catalog Number: 176379

Product name

Polyclonal

Anti-ZNF106 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Monkey

Applications

Antibody description

WB

Rabbit polyclonal antibody to ZNF106

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human zfp106. aa range:1833-1882

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:: Western Blot analysis of various cells using ZFP106 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000