

# TTC30B antibody



Catalog Number: 116379

## Product name

TTC30B antibody

## Immunogen

[Human TTC30B Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TTC30B Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse heart tissue. Observed molecular weight by Western-blot: 68-70 kDa

## Preparation

This antibody was obtained by immunization of TTC30B recombinant protein (Accession Number: NM\_152517). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116379(TTC30B Antibody) at dilution of 1:600