

# IL12A antibody



Catalog Number: 107379

## Product name

IL12A antibody

## Specificity

Human; other species not tested.

## Antibody description

IL12A Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of IL12A recombinant protein (Accession Number: NM\_000882). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107379(IL12A Antibody) at dilution of 1:2000