

BCAS1 antibody

Catalog Number: 108379

Product name

BCAS1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BCAS1 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, MCF-7 cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of BCAS1 recombinant protein (Accession Number: NM_003657). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108379(BCAS1 Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:108379(BCAS1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)