

# Collagen Type IV antibody



Catalog Number: 109379

## Product name

Collagen Type IV antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

Collagen Type IV Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human kidney tissue, rat brain tissue.

## Preparation

This antibody was obtained by immunization of Collagen Type IV recombinant protein (Accession Number: NM\_001845). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:109379(COL4A1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:109379(COL4A1 Antibody) at dilution of 1:50 (under 40x lens)