

CYP2W1 Polyclonal Antibody



Catalog Number: 161779

Product name

CYP2W1 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to CYP2W1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 141-490 of human CYP2W1 (NP_060251.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - CYP2W1 Polyclonal Antibody

Western blot analysis of extracts of mouse small intestine, using CYP2W1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.