

RLN1 Polyclonal Antibody



Catalog Number: 162779

Product name

RLN1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to RLN1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RLN1 (NP_008842.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - RLN1 Polyclonal Antibody

Western blot analysis of extracts of 22Rv1 cells, using RLN1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.