

NOD1 Polyclonal Antibody



Catalog Number: 161089

Product name

NOD1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NOD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human NOD1 (NP_006083.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - NOD1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NOD1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 120s.




Immunohistochemistry - NOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using NOD1 antibody at dilution of 1:100 (40x lens).

NOD1 Polyclonal Antibody




Catalog Number: 161089




Immunohistochemistry - NOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using NOD1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NOD1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat ovary using NOD1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - NOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using NOD1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - NOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using NOD1 antibody at dilution of 1:100 (40x lens).