

Anti-VWF antibody

Catalog Number: 102089

Product name

Polyclonal

Anti-VWF antibody

Ig Type

Immunogen

Rabbit IgG

[Human VWF \(His Tag\) recombinant protein](#)

Applications

ELISA, IHC-P

Specificity

Human vWF / von Willebrand factor

Dilutions

ELISA: 0.5-1 µg/mL

Antibody description

Rabbit polyclonal to VWF

This antibody can be used at 0.5-1 µg/mL with the appropriate secondary reagents to detect Human vWF. The detection limit for Human vWF is 0.0049 ng/well.

Preparation

Produced in rabbits immunized with purified, recombinant Human von Willebrand factor (rh vWF; Met 1-Lys 2813; NP_000543.2). vWF specific IgG was purified by human vWF affinity chromatography.

IHC-P: 0.1-2 µg/mL

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Validations


Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality



vWF Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human vWF in human tonsil with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to blood vessel.