

SEC61G antibody



Catalog Number: 115089

Product name

SEC61G antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SEC61G Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 7.7 kDa

Preparation

This antibody was obtained by immunization of SEC61G recombinant protein (Accession Number: NM_014302). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115089(SEC61G antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115089(SEC61G antibody) at dilution of 1:50 (under 10x lens)