

# Anti-DYKDDDDK-tag Mouse mAb(HRP-Conjugate)



Catalog Number: 169089

## Product name

Anti-DYKDDDDK-tag Mouse mAb(HRP-Conjugate)

## Ig Type

Mouse IgG2b

## Specificity

Tag

## Applications

WB, IP

## Antibody description

DYKDDDDK-tag Mouse mAb(HRP-Conjugate)

## Dilutions

WB: 1/1000

## Preparation

Antigen: Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (DYKDDDDK), coupled to KLH.

## Formulation

In PBS(pH 7.4) containing with 1mg/ml BSA 0.03% Proclin300 and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot analysis of extracts from CHO-K1 cells, untransfected (A) or transfected with different FLAG-fusion proteins (B,C,D), using DYKDDDDK-tag-HRP (HRP-Conjugate) Mouse mAb (1:1000 diluted).