

Human SAA4 (GST Tag) recombinant protein



Catalog Number: 501489

General Information

Gene Name Synonym

Constitutively expressed serum amyloid A protein

Protein Construction

A DNA sequence encoding the mature form of human SAA4 (P35542) (Glu19-Tyr130) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 93 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human SAA4 /GST chimera consists of 344 amino acids and has a predicted

molecular mass of 39.8KDa. It migrates as an approximately 36 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 0.5% Triton X-100, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human SAA4 Protein (GST Tag) SDS-PAGE