

Human DCXR / HCR2 (His Tag) recombinant protein



Catalog Number: 501889

General Information

Gene Name Synonym

Carbonyl reductase II; Dicarbonyl/L-xylulose reductase; Kidney dicarbonyl reductase; Short chain dehydrogenase/reductase family 20C member 1; Sperm surface protein P34H

Protein Construction

A DNA sequence encoding the mature form of human DCXR (Q7Z4W1) (Met1-Cys244) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human DCXR consists of 259 amino acids and predicts a molecular mass of 27.8 KDa. It migrates as an approximately 29 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 20% Glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human DCXR / HCR2 Protein (His Tag) SDS-PAGE