

JMY antibody

Catalog Number: 111889

Product name

JMY antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

JMY Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human colon cancer tissue, human brain tissue. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HeLa cells. Observed molecular weight by Western-blot: 111 kDa

Preparation

This antibody was obtained by immunization of JMY recombinant protein (Accession Number: BC130624). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111889(JMY Antibody) at dilution of 1:1000




IP Result of anti-JMY (IP:Catalog No:111889, 4ug; Detection:Catalog No:111889 1:600) with Jurkat cells lysate 2400ug.




Catalog Number: 111889

Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:111889(JMY Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog

No:111889(JMY Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:111889(JMY Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG