

# RHOD Polyclonal Antibody



Catalog Number: 164889

---

## Product name

RHOD Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to RHOD

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human RHOD (NP\_055393.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Immunofluorescence - RHOD Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RHOD antibody . Blue: DAPI for nuclear staining.