

TAF15 Polyclonal Antibody



Catalog Number: 165889

Product name

TAF15 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TAF15

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TAF15 (NP_631961.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - TAF15 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TAF15 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.



Immunohistochemistry - TAF15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using TAF15 antibody at dilution



of 1:100 (40x lens).

TAF15 Polyclonal Antibody



Catalog Number: 165889

Immunohistochemistry - TAF15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using TAF15 antibody at dilution of 1:100 (40x lens).