

# Anti-PLK1 Rabbit antibody



Catalog Number: 168889

## Product name

Anti-PLK1 Rabbit antibody

## Specificity

Human

## Antibody description

PLK1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human PLK1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Validations



Western blot detection of PLK1 in Jurkat,C6 cell lysates using PLK1 Rabbit mAb(1:1000 diluted).Predicted band size:68kDa.Observed band size:68kDa.