

Mouse Mmp9/Clg4b recombinant protein



Catalog Number: 518889

General Information

Gene Name Synonym

Matrix metalloproteinase-9; MMP-9; 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB

Protein Construction

A DNA sequence encoding the Mouse Mmp9/Clg4b Ala20-Pro730 is expressed with a 10His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Mmp9/Clg4b recombinant protein