

DENND3 antibody



Catalog Number: 109889

Product name

DENND3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DENND3 Rabbit Polyclonal antibody. Positive WB detected in mouse stomach tissue. Observed molecular weight by Western-blot: 136-146 kDa

Preparation

This antibody was obtained by immunization of DENND3 recombinant protein (Accession Number: BC133023). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:109889(DENND3 Antibody) at dilution of 1:1000