

NIPSNAP3A antibody



Catalog Number: 113189

Product name

NIPSNAP3A antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

NIPSNAP3A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, COLO 320 cells, HEK-293 cells. Positive IP detected in mouse brain tissue. Positive IHC detected in human colon cancer tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 28-30 kDa

Validations




Preparation

This antibody was obtained by immunization of NIPSNAP3A recombinant protein (Accession Number: NM_015469). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113189(NIPSNAP3A antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113189(NIPSNAP3A antibody) at dilution of 1:100 (under 10x lens)

Applications



ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

NIPSNAP3A antibody

Catalog Number: 113189



IP Result of anti-NIPSNAP3A (IP:Catalog No:113189, 4ug; Detection:Catalog No:113189 1:1000) with mouse brain tissue lysate 4000ug.

Immunofluorescent analysis of HEK-293 cells using Catalog No:113189(NIPSNAP3A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)