

SENP5 Polyclonal Antibody



Catalog Number: 164189

Product name

SENP5 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to SENP5

Preparation

Antigen: A synthetic Peptide of human SENP5

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SENP5 Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using SENP5 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.