

SMARCD1 antibody



Catalog Number: 107589

Product name

SMARCD1 antibody

Specificity

Human; other species not tested.

Antibody description

SMARCD1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of SMARCD1 recombinant protein (Accession Number: BC009368). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107589(SMARCD1 antibody) at dilution of 1:1000