

# MCL1L-specific antibody

Catalog Number: 106989

## Product name

MCL1L-specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MCL1L-specific Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human brain tissue, human breast cancer tissue, human cervical cancer tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human testis tissue, human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_182763). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:112550(Mcl-1l-specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:112550(Mcl-1l-specific Antibody) at dilution of 1:50 (under 40x lens)