

XKR6 antibody



Catalog Number: 116989

Product name

XKR6 antibody

Specificity

Human; other species not tested.

Antibody description

XKR6 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, Y79 cells. Observed molecular weight by Western-blot: 41-45 kDa

Preparation

This antibody was obtained by immunization of XKR6 recombinant protein (Accession Number: BC024146). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116989(XKR6 Antibody) at dilution of 1:1000