

# AMD1 antibody

Catalog Number: 107989

---

## Product name

AMD1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

AMD1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer. Observed molecular weight by Western-blot: 38kd or 32kd

## Preparation

This antibody was obtained by immunization of AMD1 recombinant protein (Accession Number: NM\_001634). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:107989(AMD1 antibody) at dilution of 1:50 (under 10x lens)