

# LNPEP antibody

Catalog Number: 112289

## Product name

LNPEP antibody

## Specificity

Human; other species not tested.

## Antibody description

LNPEP Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, COLO 320 cells, human heart tissue. Positive IHC detected in human heart tissue, human kidney tissue. Observed molecular weight by Western-blot: 117-125 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005575). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112289(LNPEP antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human heart using Catalog No:112289(LNPEP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:112289(LNPEP antibody) at dilution of 1:50 (under 40x lens)