

# Anti-JNK2 Rabbit antibody



Catalog Number: 167289

## Product name

Anti-JNK2 Rabbit antibody

WB: 1/1000-1/2000

IHC: 1/20-1/50

## Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

## Antibody description

JNK2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human JNK2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, FC, IP

## Dilutions

## Validations



Western blot detection of JNK2 in C6,3T3 cell lysates using JNK2 Rabbit pAb(1:1000 diluted). Predicted band size:48kDa. Observed band size:54kDa.



Immunohistochemistry of JNK2 in paraffin-embedded Human colon cancer tissue using JNK2 Rabbit pAb at dilution 1/20