

# HRP-conjugated GAPDH antibody



Catalog Number: 107289

## Product name

HRP-conjugated GAPDH antibody

## Specificity

Human, Mouse, Rat, Zebra Fish; other species not tested.

## Antibody description

HRP-conjugated GAPDH Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM\_001289745). Purification method: Protein A Purified.

## Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

## Validations



Western blot analysis of GAPDH in various tissues and cell lines using Proteintech antibody HRP-60004 at a dilution of 1:10000.