

# FLJ25006 antibody



Catalog Number: 110689

## Product name

FLJ25006 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FLJ25006 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of FLJ25006 recombinant protein (Accession Number: BC045622). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110689(FLJ25006 Antibody) at dilution of 1:600