

USP5 antibody

Catalog Number: 116689

Product name

USP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

USP5 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, DU 145 cells, human colon tissue, human heart tissue, MDA-MB-453s cells, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 96-100 kDa

Preparation

This antibody was obtained by immunization of USP5 recombinant protein (Accession Number: NM_003481). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116689(USP5 antibody) at dilution of 1:1000



IP Result of anti-USP5 (IP:Catalog No:116689, 4ug; Detection:Catalog No:116689 1:1000) with mouse brain tissue lysate 4000ug.