

# ACP2 antibody

Catalog Number: 107689

---

**Product name**

Polyclonal

ACP2 antibody

**Ig Type****Specificity**

Rabbit IgG

Human,Mouse,Rat; other species not tested.

**Applications****Antibody description**

ELISA, WB

ACP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HepG2 cells, HL-60 cells, mouse brain tissue. Observed molecular weight by Western-blot: 63 kDa,52 kDa

**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Preparation**


This antibody was obtained by immunization of ACP2 recombinant protein (Accession Number: NM\_001610). Purification method: Antigen affinity purified.

**Validations****Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107689(ACP2 antibody) at dilution of 1:500