

# C1R antibody

Catalog Number: 108689

---

## Product name

C1R antibody

## Specificity

Human; other species not tested.

## Antibody description

C1R Rabbit Polyclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 85 kDa, 30 kDa

## Preparation

This antibody was obtained by immunization of C1R recombinant protein (Accession Number: NM\_001733). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:108689(C1R Antibody) at dilution of 1:1000