

# HDHD3 antibody



Catalog Number: 111389

---

## Product name

HDHD3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HDHD3 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of HDHD3 recombinant protein (Accession Number: NM\_001304510). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111389(HDHD3 antibody) at dilution of 1:500