

NSUN5 antibody



Catalog Number: 113389

Product name

NSUN5 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NSUN5 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human kidney tissue, mouse skin tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human skin cancer tissue. Observed molecular weight by Western-blot: 47kd

Preparation

This antibody was obtained by immunization of NSUN5 recombinant protein (Accession Number: NM_148956). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113389(NSUN5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human skin cancer using Catalog No:113389(NSUN5 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded

NSUN5 antibody

Catalog Number: 113389

human skin cancer using Catalog
No:113389(NSUN5 antibody) at dilution of 1:50
(under 40x lens)



Immunofluorescent analysis of HepG2 cells, using
NSUN5 antibody Catalog No:113389 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).