

Syntaphilin antibody



Catalog Number: 115789

Product name

Syntaphilin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Syntaphilin Rabbit Polyclonal antibody. Positive WB detected in mouse trachea tissue. Observed molecular weight by Western-blot: 75kd

Preparation

This antibody was obtained by immunization of Syntaphilin recombinant protein (Accession Number: XM_047440615). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse trachea tissue were subjected to SDS PAGE followed by western blot with Catalog No:115789(SNPH antibody) at dilution of 1:800