

Anti-UBE3A Rabbit antibody



Catalog Number: 167789

Product name

Anti-UBE3A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

UBE3A Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human UBE3A

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of UBE3A in A549,HL-60,U2OS cell lysates using UBE3A Rabbit mAb(1:1000 diluted).Predicted band size:101kDa.Observed band size:101kDa.