

ADH6 antibody

Catalog Number: 107789

Product name

ADH6 antibody

Specificity

Human; other species not tested.

Antibody description

ADH6 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 39kd

Preparation

This antibody was obtained by immunization of ADH6 recombinant protein (Accession Number: BC039065). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107789(ADH6 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using ADH6 antibody Catalog No:107789 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).