

Anti-MRPS21 antibody



Catalog Number: 176099

Product name

Anti-MRPS21 antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Rabbit polyclonal antibody to MRPS21

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrps21. aa range:38-87

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.;; Western Blot analysis of various cells using MRP-S21 Polyclonal Antibody diluted at 1: 500