

# Anti-HA-Tag Mouse mAb (HRP Conjugate)



Catalog Number: 169099

## Product name

Mouse IgG1

Anti-HA-Tag Mouse mAb (HRP Conjugate)

## Applications

## Antibody description

WB, ELISA

HA-Tag Mouse mAb (HRP Conjugate)

## Dilutions

## Preparation

WB: 1/500-1/5000

Antigen: YPYDVPDYA conjugated to KLH.

## Validations

## Formulation

In PBS(pH 7.4) containing with 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot analysis of extracts from CHO-K1 cells (B) or CHO-K1 cells transfected with HA-fusion proteins (A), using HA-Tag Mouse mAb (HRP Conjugate) (1:1000 diluted).