

Anti-HA-Tag Mouse mAb (HRP Conjugate)



Catalog Number: 169099

Product name

Mouse IgG1

Anti-HA-Tag Mouse mAb (HRP Conjugate)

Applications

Antibody description

WB, ELISA

HA-Tag Mouse mAb (HRP Conjugate)

Dilutions

Preparation

WB: 1/500-1/5000

Antigen: YPYDVPDYA conjugated to KLH.

Validations

Formulation

In PBS(pH 7.4) containing with 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot analysis of extracts from CHO-K1 cells (B) or CHO-K1 cells transfected with HA-fusion proteins (A), using HA-Tag Mouse mAb (HRP Conjugate) (1:1000 diluted).