

EXO1 antibody

Catalog Number: 110499

Product name

EXO1 antibody

Specificity

Human; other species not tested.

Antibody description

EXO1 Rabbit Polyclonal antibody. Positive IP detected in mouse thymus tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 115-120 kDa

Preparation

This antibody was obtained by immunization of EXO1 recombinant protein (Accession Number: XM_011544322). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110499(EXO1 antibody) at dilution of 1:800



IP Result of anti-EXO1 (IP:Catalog No:110499, 4ug; Detection:Catalog No:110499 1:500) with mouse thymus tissue lysate 2800ug.