

# PIGZ antibody

Catalog Number: 113899

## Product name

PIGZ antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PIGZ Rabbit Polyclonal antibody. Positive IHC detected in human colon. Positive WB detected in human brain tissue, HEK-293 cells, human testis tissue. Observed molecular weight by Western-blot: 63kd

## Preparation

This antibody was obtained by immunization of PIGZ recombinant protein (Accession Number: NM\_025163). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113899(PIGZ antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon using Catalog No:113899(PIGZ antibody) at dilution of 1:100 (under 10x lens)