

# YWHAZ antibody



Catalog Number: 116899

## Product name

YWHAZ antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

YWHAZ Rabbit Polyclonal antibody. Positive WB detected in Raji cells, Jurkat cells, mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 28-30 kDa

## Preparation

This antibody was obtained by immunization of YWHAZ recombinant protein (Accession Number: XM\_005251061). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116899(YWHAZ Antibody) at dilution of 1:1000