

# Anti-PGM1 Rabbit antibody

Catalog Number: 168899

## Product name

Anti-PGM1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PGM1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human PGM1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


1:2000

1:100


1:20

1:20

## Validations



Western blot detection of PGM1 in C6,3T3,HeLa cell lysates using PGM1 Rabbit mAb(1:1000 diluted).Predicted band size:62kDa.Observed band size:62kDa.



Immunocytochemistry of PGM1(green) in HeLa cells using PGM1 Rabbit mAb at dilution 1/50, and DAPI(blue)