

HUS1B Polyclonal Antibody



Catalog Number: 162199

Product name

HUS1B Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to HUS1B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human HUS1B (NP_683762.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Immunohistochemistry - HUS1B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using HUS1B antibody at dilution of 1:100 (40x lens).