

# CAPN6 Monoclonal Antibody



Catalog Number: 164199

## Product name

CAPN6 Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to CAPN6

## Preparation

Antigen: A synthetic peptide of human CAPN6

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - CAPN6 Monoclonal Antibody

Western blot analysis of extracts of human placenta, using CAPN6 antibody .Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.