

# PRKCDBP antibody

Catalog Number: 114199

## Product name

PRKCDBP antibody

WB: 1:200-1:1000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

PRKCDBP Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive IF detected in MCF-7 cells. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, MCF7 cells. Observed molecular weight by Western-blot: 35-40kd

## Validations



## Preparation

This antibody was obtained by immunization of PRKCDBP recombinant protein (Accession Number: NM\_145040). Purification method: Antigen affinity purified.

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114199(PRKCDBP antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IF, IP

## Dilutions


Recommended Dilution:



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114199(PRKCDBP antibody) at dilution of 1:100 (under 40x lens)

# PRKCDBP antibody

Catalog Number: 114199



Immunofluorescent analysis of MCF-7 cells, using PRKCDBP antibody Catalog No:114199 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IgG (red).

IP Result of anti-PRKCDBP (IP:Catalog No:114199, 4ug; Detection:Catalog No:114199 1:600) with mouse lung tissue lysate 2800ug.

